



Website: www.somamedical.net


Sales Inquiry

AccuRead Arm DB61M (Blood Pressure Monitor with WHO Light Sign Indicator)



Feature:

- Fully automatic inflation and deflation
- Super large LCD display (size: 6 x 7 cm)
- 4 person use/99 memory each
- Irregular heartbeat detector
- 4 alarm clock settings
- 5 mode average value of reading
- 6 colour WHO (World Health Organization) light indicator



**WHO (World Health Organization)
light sign indicator**

Specifications

Blood Pressure

Measurement Range: 20~280mm Hg
Accuracy: +/- 3mmHg
Display: Systole, Diastole, Pulse

Pulse

Measurement Range: 40~195 pulse/min.
Pulse Rate Accuracy: +/- 5% of the reading

Power Dimension

Power: 3V, 4 x AAA size
Dimension: 4.8" (W) x 6" (H) x 5.3" (D)
Weight: Approx. 344g (excl. batteries), cuff 120g

Arm Circumference : 22 - 34 cm (Special order for extra size)

Auto Power off after 3 minutes

Product Qualification

CE0434: This product fulfils the requirement for Directive 93/42/EEC concerning Medical Device.

FDA KO14141/ISO 9001/MDD Annex V: Manufactured in accordance with International standards

TRULY